



PORTABLE DCM CABINET

DESCRIPTION

Portable cabinet for wall installation, hanging on the wall and manual transport, with a maximum capacity to house two or three DCM100 and / or DCM200 medical devices. It has an interchangeable base for two or three units.

INFORMATION

SKU	CATEGORY	WEIGHT	DIMENSIONS	MATERIAL	COLOR
ACC3	ACCESSORIES DMC	500 grs.	L: 26 cm W: 8 cm H: 21 cm	ECO-LEATHER POLYURETHANE	GREEN



- Shock resistance.
- Weight support: Up to 20 Kg.
- Wash exclusively by hand, with a damp cloth and at room temperature. Top magnetic closure.
- 1 Fixed strap for manual transport.
- It has two holes in the back for fixed installation to the wall and a hook in the upper part for semi-fixed installation to the wall.