



DCM100 - NEBULIZER 100 ml

DCM200 - NEBULIZER 200 ml

## DCM

### HAND AND FACE DECONTAMINATOR

ONLY USE IN CASE OF CHEMICAL BURN

#### DESCRIPTION

Face and Hand Decontaminant with Anphoterol® solution, it is a class II medical device, portable and easy to apply. It is an aerosol washing system with first aid fogging action, for accidents caused by splashing with chemical agents of acid or alkaline pH. **Presented in two formats of 100 ml and 200 ml of Anphoterol®** solution, indicated for affected areas of up to 3% and 9% of the total body surface respectively. Anphoterol® solution is a non-irritating, non-allergic, biodegradable, non-toxic product (to the skin and eyes) since its compounds are present in the human body.



#### INFORMATION

SKU	CONTAIN	DIMENSIONS	PROPELLANT GAS	EXPIRATION DATE	USE	WASH TIME
DCM100	100 ml 104 g	Ø 53 mm H: 180 mm	Bag-on-Valve system	24 months	HANDS / EYES FOOT / FACE	70 sec.
DCM200	200 ml 208 g	Ø 53 mm H: 210 mm			ARMS / HANDS FEET / LEG EYES / HEAD	140 sec.

- **Conforms to ISP DDM / 2027/19.**
- Simple wash control: complete and single use of the device.
- Mist and moisturizing spray.
- Class II medical device.
- With a biodegradable solution developed with compounds present in the human body.
- Uniform washing in the affected area.
- It does not cause damage or trauma to the skin or eye area.
- Light and easy to handle presentation.
- Indicated for affected areas of up to 3% (DCM100) and 9% (DCM200) of the total body surface.

