

EMERGENCY CASE

DESCRIPTION

First aid kit for chemical and thermal burns. This complete Kit gathers elements to attend the first aid of chemical and temperature accidents; in a practical anti-impact, waterproof, black case with white reflective tape around. Inside it has transparent pvc cases with closure, which allow a previous distribution of the elements according to the type of emergency that is required to attend. On the back you will find an external pocket for storing documents.



INFORMATION

SKU	WEIGHT	DIMENSIONS	MATERIAL	COLOR
KIT1	2,35 Kg.	H:27 cm- L:35 cm W:10 cm	BODY: PVC STRAP: NYLON	BLACK

- Internal padded anti-impact reinforcement to protect the contents against possible falls of the briefcase.
- Case covered with transparent PVC easy to clean. Back zipper.
- Reflective tape for immediate viewing.

CONTAIN

CHEMICAL BURNS

3 DCM 100 nebulizers.

ANPHOTEROL

1 DCM 200 nebulizer.

BURN TREATMENT

- 3 Hydrogel dressings for burns 10x10cm.
- 1 Hydrogel dressing for burns 20x20cm.
- 1 Hydrogel dressing for burns 30x40 cm Face.
- 15 Hydrogel sachet 3.5 gr.
- 1 Hydrogel spray 120 ml.

BASIC KIT - FIRST AID

- 1 Extrication scissors tip duck 15 cm
- 2 Elasticated gauze bandages 8cm x 4mts.
- 5 Adhesive tape
- 1 Pair of vinyl gloves for procedures.
- 1 Emergency briefcase and documentation.