



DMT2L DECONTAMINATOR - 2 lts

DMT5L DECONTAMINATOR - 5 lts

## DMT


### MEMBER DECONTAMINATOR AND TORSO

ONLY USE IN CASE OF CHEMICAL BURN

### DESCRIPTION

Limb and Torso Decontaminant with Anphoterol® solution, it is a class II medical device, portable and easy to apply. It is a portable washing shower with nebulizing action presented in **three formats of 2 liters and 5 liters of Anphoterol®** solution, indicated as a first aid system for accidents caused by splashes with chemical agents of acid or alkaline pH in large areas of the body, such as head, arms, legs or torso.

Anphoterol® solution is a non-irritating, non-allergic, biodegradable, non-toxic product (to the skin and eyes) since its compounds are present in the human body.



### INFORMATION

SKU	CONTAIN	DIMENSIONS	PROPELLANT GAS	EXPIRATION DATE	USE	WASH TIME
DMT2L	2 L 2,1 kg	Ø 100 mm H: 450 mm	14 kg / cm <sup>2</sup> Nitrogen	24 months	ARM LEG TORSO	-4 min.
DMT5L	5 L 5,2 kg	Ø 115 mm H: 520 mm				-8 min.

- **Conforms to ISP DDM / 2027/19**
- Simple wash control: Complete and single use of the device.
- Easy to differentiate, shiny and made of stainless steel.
- Transportable from one work area to another.
- Class II medical device.
- With a biodegradable solution developed with compounds present in the human body.
- Mist and moisturizing spray.
- 1.2 m long hose for easy maneuvering.
- Seal and security pin to guarantee control of use.
- Accessories included for wall installation.
- Indicated for affected areas on 9% of the surface total body.